ENGLISH

Medicel ACCUJECT™ / ACCUJECT™ SCREW
Single-use Lens Injection System for injection of ONE-PIECE FOLDABLE LENSES

DESCRIPTION
With the aid of the ACCUJECT™ (Fig. 1a) and ACCUJECT™ SCREW (Fig. 1b) injection systems, the majority of all foldable, 1-piece hydrophobic acrylic lenses

with an optical diameter up to 6.5 millimetres can be injected, provided that the safe injection of the lens using the respective ACCUJECT™ or ACCUJECT™ SCREW injection system was validated in advance by the lens manufacturer as per ISO 11979. The cartridges of the ACCUJECT™ and ACCUJECT™ SCREW models are designed with a large covered loading chamber. This enables the lenses to be loaded quickly, precisely and easily.

APPLICATION
Folding and injection of 1-piece folding lenses in the capsular bag or the sulcus after extracapsular cataract extraction.

INSTRUKTIONEN FÜR DIE NUTZUNG
1. Open the blister in a sterile environment and remove the sterile ACCUJECT™ or ACCUJECT™ SCREW injector set.

2. Fill the cartridge nozzle with viscoelastic solution of low to moderate viscosity directly from the cartridge tip (Fig. 2). Now wet the bottom of the loading chamber with viscoelastic solution from the end of the loading chamber cover to make this slippery (Fig. 3). Do not completely fill the chamber with viscoelastic solution, as this can move the lens during insertion. Also put one drop of viscoelastic solution on the silicone tip. When using certain hydrophobic lens materials, the cartridge tip and loading chamber should be completely rinsed out with a BSS afterwards in order bag to achieve maximum hydration. This step is generally not necessary when using hydrophilic lens materials. The exclusive use of a BSS solution, or the use of viscoelastic solutions, with high viscosity is explicitly not recommended. The lens should be injected within 3 minutes after

filling viscoelastic solution. Viscoelastic substances can lose their lubrication properties if they are exposed to air for too long.

3. Put the lens down on the rear, uncovered part of the loading chamber. Use the angled forceps to press against the rear end of the optic.

IMPORTANT: In the case of lenses with c-loop haptics, press only against the optic and not on the rear haptic.

In the case of a loading chamber with cover as shown in Fig. 4a, push the lens into the loading chamber until the forceps butt against the end of the guide slot of the top cover (Fig. 4b). Hold the forceps in this position and use the other hand to close the wings of the loading chamber (Fig. 4c).

IMPORTANT: When using lenses with c-loop haptics, pull the forceps out only after the wings of the loading chamber

have already been somewhat pressed together and pressure has already been exerted on the lens. This prevents the lens from slipping back out of its pushed-in position and the pre-folded front haptic from stretching back to its original shape.

In the case of a loading chamber with cover as shown in Fig. 5a, slide the lens into the loading chamber so that the rear haptic is completely underneath the cover. Push the injector plunger forwards until the lateral marking notch of the plunger is flush with the end of the injector housing and the plunger engages with a light "click" in the injector housing (Fig. 5b).

4. Close the cartridge by pressing the wings together (Fig. 6). As soon as the "click/lock" mechanism clicks in place, the lens is safely loaded and ready for injection.

IMPORTANT: Pressing the wings together prematurely (even

just slightly) **before sliding the lens in** can cause damage to the system and affect the functioning of the injector.

5. Press the injector plunger forward until the blue silicone tip of the ACCUJECT™ plunger has reached the end of the injector housing (Fig. 7), or until the threads of the ACCUJECT™ SCREW injector are tangibly touching the injector housing.

IMPORTANT: Pull the plunger back a few millimetres and then push it forward once again. This step ensures that the rear haptic is not jammed.

6. Guide the point of the cartridge through the incision and push it over the iris to the near edge of the pupil.

7. Slowly push the ACCUJECT™ injector plunger forward, or with the ACCUJECT™ SCREW injector, slowly turn it in a clockwise direction to push the lens forward.

8. Depending on the lens model, support the correct ejection of the lens by slightly turning the injector. Slowly inject the lens into the eye and simultaneously withdraw the instrument from the eye. To avoid swelling of the silicone tip during ejection from the cartridge, only push the plunger until the lens has emerged completely, even if the plunger is not yet at the limit.

9. If the rear haptic is jammed between the silicone tip and the cartridge tip after the plunger has been completely pushed in, retract the plunger together with the silicone tip a few millimetres until the rear haptic is exposed; then inject the rear haptic in a second thrust.

10. Carefully position the lens, rotating it if necessary with the aid of using a suitable positioning hook.

11. Remove the viscoelastic material thoroughly from the eye

and the lens with the standard irrigation and aspiration techniques.


PLEASE ENSURE: Do not re-sterilize / prepare ACCUJECT™ or ACCUJECT™ SCREW injection System. Single use only. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

OPERATIONAL PROCEDURE
The appropriate surgical techniques are the responsibility of the respective surgeon. He or she must assess the appropriateness of the relevant procedure based on his or her education and experience.

GUARANTEE AND LIABILITY LIMITATION

The manufacturer guarantees that this product has been manufactured with the appropriate care and assumes no liability for side effects or resulting damages, losses or costs that may arise as a result of the direct or indirect use of this product. Manufacturer's liability is restricted to the performance of repairs resulting from product defects, which are clearly not the result of incorrect handling or the use of lenses not validated with this injector model.

ATTENTION: US federal regulations restrict the sale of this product to medical practitioners and those acting on their behalf.*

ACCUJECT™
ACCUJECT™ SCREW
— backloaded —

medicel
SWISS TECHNOLOGY FOR SURGERY

Instructions for use
Gebrauchsanweisung
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso

MADE IN SWITZERLAND

Manufacturer:
MEDICEL AG
Dornierstrasse 11
CH-9423 Altenrhein
SWITZERLAND

Tel. +41 71 727 10 50
Fax +41 71 727 10 55
info@medicel.com
www.medicel.com

medicel
SWISS TECHNOLOGY FOR SURGERY

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






ENGLISH

Medicel
ACCUJECT™ SCREW
Single-use Lens Injection System
 for injection of
ONE-PIECE FOLDABLE LENSES

DESCRIPTION
 With the aid of the ACCUJECT™ SCREW (Fig. 1) injection systems, the majority of all foldable, 1-piece acrylic lenses with an optical diameter up to 6.5 millimetres can be

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medicel
 SWISS TECHNOLOGY FOR SURGERY

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MADE IN SWITZERLAND

Manufacturer:
MEDICEL AG
 Dornierstrasse 11
 CH-9423 Altenrhein
 SWITZERLAND

Tel. +41 71 727 10 50
 Fax +41 71 727 10 55
 info@medicel.com
 www.medicel.com

injected, provided that the safe injection of the lens using the ACCUJECT™ SCREW injection system was validated in advance by the lens manufacturer as per ISO 11979. In contrast to traditional injection systems, the cartridge of the ACCUJECT™ system is integrated into the injector. The cartridges of the ACCUJECT™ SCREW models are designed with a large loading chamber. This enables the lenses to be loaded quickly, precisely and easily from above.

APPLICATION
 Folding and injection of 1-piece folding lenses in the capsular bag or the sulcus after extracapsular cataract extraction.

INSTRUCTIONS FOR USE

- Open the blister in a sterile environment and remove the sterile ACCUJECT™ SCREW injector set.
- Fill the cartridge nozzle with viscoelastic solution of low to

moderate viscosity directly from the cartridge tip. (Fig. 2) Fill the ACCUJECT™ SCREW loading chamber sufficiently with viscoelastic solution (Fig. 3). In addition apply a drop of the viscoelastic solution directly to the silicone tip. When using certain hydrophobic lens materials, the cartridge tip and loading chamber should be completely rinsed out with a BSS afterwards in order to achieve maximum hydration. This step is generally not necessary when using hydrophilic lens materials. The exclusive use of a BSS solution, or the use of viscoelastic solutions, with high viscosity is explicitly not recommended. The lens should be injected within 3 minutes after filling viscoelastic solution. Viscoelastic substances can lose their lubrication properties if they are exposed to air for too long.

- Hold the ACCUJECT™ SCREW injection system in

such a manner that the rear cartridge wing can be guided with the index finger of your left hand. Open the loading chamber sufficiently wide so that the lens can be inserted easily (Fig. 4).

- Position the lens in the middle of the loading chamber between the silicone tip and cartridge wall (Fig. 5).
- Close the cartridge and exert slight pressure onto the lens optic, using a forceps (Fig. 6). This will ensure that the lens is folded in the right direction. Once the "click-lock" mechanism engages the lens is securely loaded and ready for injection (Fig. 7).
- Press the injector plunger forward until the threads of the injector are tangibly touching the injector housing.
- Guide the point of the cartridge through the incision and push it over the iris to the near edge of the pupil.

- Slowly turn the injector plunger in a clockwise direction to push the lens forward.
- Depending on the lens model, support the correct ejection of the lens by slightly turning the injector. Slowly inject the lens into the eye and simultaneously withdraw the instrument from the eye. To avoid swelling of the silicone tip during ejection from the cartridge, only turn the plunger until the lens has emerged completely, even if the plunger is not yet at the limit.
- If the rear haptic is jammed between the silicone tip and the cartridge tip after the plunger has been completely pushed in, retract the plunger together with the silicone tip a few millimetres until the rear haptic is exposed; then inject the rear haptic in a second thrust.
- Carefully position the lens, rotating it if necessary with the

aid of using a suitable positioning hook.

- Remove the viscoelastic material thoroughly from the eye and the lens with the standard irrigation and aspiration techniques.

PLEASE ENSURE: Do not re-sterilize/prepare the instrument from the direct or indirect use of this product. Manufacturer's liability is restricted to the performance of repairs resulting from product defects, which are clearly not the result of incorrect handling or the use of lenses not validated with this injector model.

OPERATIONAL PROCEDURE
 The appropriate surgical techniques are the responsibility of the respective surgeon. He or she must assess the appropriateness of the relevant procedure based on his or her education and experience.

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