ENGLISH

Medichel VISCOJECT™ MULTIJECT™

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DESCRIPTION
The VISCOJECT™ injection system permits, thanks to the VISCOJECT™ hydraulic system, the implementation of foldable intraocular lenses through incisions of about 1.5 to 2.8 mm, depending on the VISCOJECT™ model and operating technique used. With the VISCOJECT™ injection system, one-piece foldable intraocular lenses are not suited for implantation with an optical diaphragm or foldable, one-piece acrylic lenses with an optical diaphragm can be injected, provided that the provider of up to 6.5 millimeters can be injected, with the particular VISCOJECT™ injection system.

TREATMENT ACCORDING TO ISO EN 17664
General principles: Before every use, all instruments must be cleaned, disinfected and sterilized; this applies particularly to the first use after delivery since all instruments are delivered sterilized (cleaning and disinfection after removal of protective shipping packaging; sterilization). An effective cleaning and disinfection is an absolute necessity for effective sterilization. As part of your responsibility for the sterility of instruments during use, please ensure that only appropriate devices and product-specific procedures for cleaning

reuse or resterilization may reduce the product's performance, which can result in serious impairments to the patient's health and safety or even damage to said instrument.

APPLICATION
Folding and injecting of one-piece foldable lenses into the capsular bag or sacculus after extra-capsular cataract extraction. Three-piece foldable intraocular lenses are not suited for implantation with an optical diaphragm or foldable, one-piece acrylic lenses with an optical diaphragm can be injected, provided that the provider of up to 6.5 millimeters can be injected, with the particular VISCOJECT™ injection system.

CLEANING AND DISINFECTION
Principles according to ISO EN 15880: For the cleaning and disinfection of instruments, a mechanical process (disinfection) should be used, if possible. Due to its clearly reduced effectiveness, a manual procedure – even when using an ultrasound bath – should only be used if the disinfectant used is for protecting persons and must not replace the pretreatment step to be carried out later, after cleaning has been done. Check all instruments for visible contamination. Repeat the cleaning and disinfection device or other measures must ensure sufficient drying.

INSPECTION AND TESTING
Check all instruments after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chipping and contamination and silt out damaged instruments. Instruments that are still dirty must be cleaned again.

MACHINE CLEANING/DESINFECTATION
In selecting the disinfectant, make sure that the disinfectant basically has a tested effectiveness (e.g. DGHM, or FDA approval or CE marking), a tested program (disinfection) and product-specific validation (at least 10 minutes). If the disinfectant is not chemically disinfectant (at least 10 minutes), it does not contain the following ingredients: strong acids or bases, -alcohols, -ketones, -halogens, -halogenated hydrocarbons, -furan, -phenols, -methylene chloride, -nitrobenzene.

We recommend sterilizing instruments in the sterilization tray and sterilization containers designed for this purpose; single-use packages (single or double packaging) and other sterilization containers can be used, however, as long as these meet the following requirements: compliance with DIN EN ISO/ANSI/AAMI ISO 11607 part 2+ suitable for steam sterilization (temperature resistance to at least 141 °C, sufficient vapor permeability, sufficient protection of the instruments or sterilization packagings from mechanical damages, regular maintenance according to the manufacturer guidelines (see chapter "Material Durability"). Sterilization temperature of 134 °C/+19 °C (274 °F/+62.4 °F), with saturated steam pressure of 3.41 bar, sterilization time (exposure time at sterilization temperature) of at least 20 minutes at 121 °C (250 °F) or 5 minutes at 134 °C (274 °F).

REINJECTION AND STERILIZATION
For the sterilization, only the following sterilization procedures may be used: no use of hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, and no plasma sterilization.

STORAGE
After sterilization, the instruments must be stored dry.

PLEASE NOTE:
The injector's cartridge may not be treated or resterilized. The product to medical practitioners and those acting on their behalf.

cartridge may only be used once. Information concerning the use of the cartridge is contained in the instructions for use of the VISCOJECT™ BIO and VISCOJECT™ cartridge set.

DISINFECTION
The appropriate surgical technique is the responsibility of the respective surgeon. Based on his education and experience, he is to evaluate the suitability of the respective procedure.

GUARANTEE AND LIMITATION OF LIABILITY
The manufacturer guarantees that his product was produced with appropriate care and shall assume no responsibility for incidental or consequential damages, losses or costs that should result directly or indirectly from the use of this product. Liability is solely limited to claim-related repairs that must be performed on the product which is not clearly attributable to incorrect handling or use of the product.

RELIABILITY
The instruments can be reused as long as they are handled carefully and are undamaged and uncontaminated. The user takes responsibility for the use of damaged or contaminated instruments.

STERILIZATION
For the sterilization, only the following sterilization procedures may be used: no use of hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, and no plasma sterilization.

STERILIZATION
After sterilization, the instruments must be stored dry.

MAINTENANCE
Information concerning the cleaning agents and disinfectants, make sure that the cartridge is not cleaned (the following ingredients: strong acids or bases, -alcohols, -ketones, -halogens, -halogenated hydrocarbons, -furan, -phenols, -methylene chloride, -nitrobenzene).

MAINTENANCE
It is the responsibility of the user to clean and disinfect the instruments in the sterilization tray and sterilization containers.

DISINFECTANT
We recommend sterilizing instruments in the sterilization tray and sterilization containers designed for this purpose; single-use packages (single or double packaging) and other sterilization containers can be used, however, as long as these meet the following requirements: compliance with DIN EN ISO/ANSI/AAMI ISO 11607 part 2+ suitable for steam sterilization (temperature resistance to at least 141 °C, sufficient vapor permeability, sufficient protection of the instruments or sterilization packagings from mechanical damages, regular maintenance according to the manufacturer guidelines (see chapter "Material Durability"). Sterilization temperature of 134 °C/+19 °C (274 °F/+62.4 °F), with saturated steam pressure of 3.41 bar, sterilization time (exposure time at sterilization temperature) of at least 20 minutes at 121 °C (250 °F) or 5 minutes at 134 °C (274 °F).

ATTENTION: US federal regulations restrict the sale of this product to medical practitioners and those acting on their behalf.

VISCOJECT™
MULTIJECT™

Instructions for use
Gebrauchsweisung
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso
Kullanım Talimatı

medichel
SWISS TECHNOLOGY FOR SURGERY

Instructions for use
Gebrauchsweisung
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso
Kullanım Talimatı

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