

20 PROGRAMS

Two-channel electrotherapy device for home use for the treatment of pain and the most common pathologies, in particular back pain and cervical pain. It also includes programs for beauty and sport and operates on rechargeable batteries.



16 TENS and PAIN

- TENS
- CERVICAL PAIN
- LUMBAGO
- SCIATICA
- ARTHRITIS
- OSTEOARTHRITIS
- CARPAL TUNNEL
- MUSCLE INJURY
- EPICONDYLITIS
- MUSCLE PAIN
- MENSTRUAL PAIN

4 SPORT, FITNESS, BEAUTY

- ABDOMEN MAXIMUM STRENGTH
- ABDOMEN SCULPTING



Watch our videos on [YouTube](#)

TECHNICAL FEATURES

Channels: 2 independent channels (4 electrodes) - **Frequency:** 0-150 Hz

Ampiezza impulso: 50-400 μ s - **Power:** 100 mA peak/peak

Power supply: rechargeable battery - **Size:** 8x13x2.5 cm

EQUIPMENT

- 1 carrying case
- 1 DUO TENS stimulator
- 2 cables for electrode connection
- 4 self-adhesive square electrodes
- 4 self-adhesive rectangular electrodes
- 1 charger
- 1 operating manual

