



PB 331 / PB 337 LOW FRONT LIFTING CARE BED:

Timeless front lifting bed with well-proven technology



With its beech wooden parts, the closed wooden fronts and the silver color powdered metal pieces this high class care bed blends in with its room surroundings. To avoid any sort of squeezing hazard there are side rail end caps at the ends of the side rails that minimize the space between side rails and wooden fronts in a way that prevent accidental squeezing. The side bars can be lowered and slid directly onto each other so that they do not pose any hindrance for entering or exiting the bed.

This bed complies to the highest standards. Its 4-sectional lying surface is individually adjustable with 2 independent motors. Furthermore the adjustment of the foot part is assisted by 2 ratchets. For the standard version of the bed, the lying area with wooden sprung slats and durability adjuster comes with a surface area of 90 cm x 200 cm. It can also be ordered at a width of 80 cm and in different lengths from 180 cm to a length of 220 cm. Against surcharge an easy to clean metal mesh lying area can be ordered.

The products will be delivered using transport-racks (at a maximum of 2 beds per palette). The beds have been partly preassembled and can be easily setup by a single person. The 4-sectional lying surface is electrically adjustable between 40 cm and 80 cm (upper end of the lying area). All electronic components have been manufactured according to latest research and are in accordance with all valid safety requirements and standards. Every single adjustment option of the bed is comfortably controllable by using the splash-proof manual control switch and any single function can be individually locked.

General Information:

Save work load:	175 kg
Max. User weight:	135 kg
Total weight:	102,5 kg
Heaviest single part:	Front komplett 24,9 kg
Dimensions:	Length: 215 cm, Width: 102 cm
Back part length:	80 cm
Heading adjustability:	70°
Foot adjustability:	45°
Height under mattress support:	22 cm (PB337 – 5 cm)
Caster size:	100 mm always with brake.

PB. 331 / PB 337 full fills all valid safety requirements and standards according to EU guidelines 93/42/EWG for devices, complies medical with the CE mark of conformity.

loading