

**New**

# EL701 LED Examination Light

- \*Perfect for the ENT/GYN/Dermatology/Dental office or for specialties .
- \*With a high performance LED, Provides you with an absolutely homogeneous field of light .
- \*Touchless on/off
- \* With a High-capacity rechargeable battery , could be used when it is charging.
- \*1.2 m goose neck is as long as conventional pipes, Two-segment design could be bent radius to any angle and in all directions.



EL701

## EL701 series

Optical specificatoinis	
Diameter of the light spot at a distance of 300mm	30-100mm
Illumination at a distance of 300mm	20000lux
Colour temperature	≥4000K
Colour rendering index Ra	≥90
Life time of the light	50,000 hours
Misc data	
Mounting options	Wheeled stand
	Wall mount
	Table-top mounting
Expected service life	≥5 years

Examination Light	
Order Code	Description
<b>701E2000</b>	EL7010 LED examination light with clamp for table-top mount
<b>701E3000</b>	EL7010 LED examination light with wall mount
<b>701E5000</b>	EL7010 LED examination light with wheeled stand
<b>701E2050</b>	EL7015 LED examination light with clamp for table-top mount,with battery
<b>701E5050</b>	EL7015 LED examination light with wheeled stand,with battery



EL701