

Maximum load:

Active Lift: 175 kg

Material: powder-coated steel

Weight: Total: 43.0 kg

Heaviest removable part: 22.9 kg

Chassis:

Lowest point: 60 mm

The highest point: 120 mm

Width closed: 575 mm

Maximum spread: 1175 mm

Overall length: 1185 mm

Maximum height in the extended arm: 2040 mm

Lowest height support arm: 635 mm

Castors:

Standard front: 75 mm double castors

Standard rear: 75 mm single castor with brake

Turning circle: 1063 mm

Emergency lowering: Mechanical and Electrical

Lifting speed (without load):

Slingbar 350: 54 mm / s

Maximum noise level: 42 dB (A)

Protection class: IP X4

