

# Princess II Water Birth Pool



The compact yet surprisingly Princess II is the perfect choice for rooms with limited space.


This perfectly formed gem of a pool now features a specially designed labour seat to support mothers in upright postures and give them a place to rest and welcome their baby.

Though we've pared things down to the minimum to create this highly utilitarian pool we made sure to include all of the essentials that make our pools stand out from the rest.



## Specifications

Length:	1600mm
Width:	1100mm
Height:	770mm
Depth:	640mm
Weight (empty):	70kg
Capacity:	425 litres
Material:	Ficore® composite
Guarantee:	Lifetime



# Princess II Water Birth Pool

ACTIVE  
BIRTHPOOLS

Multi-colour LED Lighting System gives mothers six ambient colours to choose from

The ultra-hard material is less slippery and more resistant to bacterial and cleansing agents

The key-hole shape channel in the rim protects mothers from hitting their head on the spout

Bluetooth Sound System with integral speakers for mothers to play music of their choice

The Labour Support Seat provides support in upright positions and a place for mothers to rest

Bespoke drainage system with pop-up waste and rim mounted control

Midwives sit in comfort with their legs under the pool and arms resting on the extra-wide rim

The extra-wide front rim makes it easy for mothers to get in and out of the pool

The bull-nose shape wrap-around rim provides mothers and midwives with maximum comfort and support

The Lumbar Support provides exceptional support for mothers in squatting positions

Safe, practical options to facilitate emergency evacuation with or without a portable hoist

There's plenty of room for mothers to move freely to gain maximum benefit from water immersion

