

KINETEC™ Maestra™ Hand & Wrist CPM Machine

Rehabilitation



Post-Operatively



The **CPM Solution** handling every single hand and wrist pathology

- **A useful tool for protocols**
dedicated to surgeons,
physiotherapists and occupational therapists
- **Rehabilitation/Post-operatively**
- **Each hand and wrist movement possible**
- **Analytic / Ambulatory**



Surgical techniques, users information and protocols available on CD ROM.

Ordering information "Free of Charge" : 4696718097

Indications

- Flexor and extensor tendon tenolysis
- Aponeurectomies for Dupuytren's disease
- Open reduction and internal fixation of intra-articular, diaphyseal, metaphyseal and epiphyseal fractures of the phalanges
- MCP arthrolysis
- Prosthetic replacement of the MCP, PIP and DIP joints
- Rheumatoid/neurological or post-burn stiffness
- Capsulotomies, arthrolysis and tenolysis for post-traumatic stiffness of MCP, PIP and DIP joints
- Synovectomies



REHABILITATION

Maestra
9 Movements
Hand & Wrist



- 1• Composite fist (MCP-PIP-DIP)
-30° to 225°



- 2• Ulnar deviation
-30° to 60°
- 3• Extension / Flexion of the wrist
-50° to 90°
- 4• Extension / Flexion of the wrist with MCP
-50° to 140°
- 5• Pronation / supination
-90° to 90°



- 6• Thumb Opposition
0° to 180°



- 7• Intrinsic +
0° to 90°
- 8• Intrinsic -
0° to 180°
- 9• Selective mobilisation of the DIP joint
0° to 70°

POST-OPERATIVELY

Maestra Portable
3 Movements
Longs Fingers



- 1• Composite fist (MCP-PIP-DIP)
-15° to 270°



- 2• Intrinsic -
0° to 180°



- 3• Intrinsic +
0° to 90°

Ordering Information

- 4621005502 Complete apparatus • Maestra Portable
- 4650001074 Set of 4 finger attachments



Ordering Information

- 4621005002 Complete apparatus • Maestra
- 46EXTGAR3G 3-Year Warranty Extension
- 46500001793 Complet exchange pads
- 4650001470 Exchange arm rest strap support (x2)
- 4650000753 Exchange arm rest support pads (x1)
- 4650001462 Exchange Ulnar deviation support pads (x1)

