

DIACARE 7000

TECARTHERAPY

PROFESSIONAL DEVICE TO TAKE CARE OF THE BODY USING THE INTERNAL ENERGIES OF THE BODY AND THE NATURAL AUTOGENERATION MECHANISMS OF TISSUES AND MUSCLES. THE COLOR TOUCH SCREEN NAVIGATION SYSTEM MAKES THE RESEARCH AND THE PARAMETER MODIFICATION EASIER. THE COMPLETE SET OF PROGRAMS, THE ATLAS THAT DESCRIBES WITH PICTURES THE TREATMENTS IN THE DIFFERENT PARTS OF THE BODY AND A SECTION WHERE THE OPERATOR CAN STORE THE TREATMENTS ARE JUST A FEW OF THE REASONS THAT MAKE DIACARE 7000 THE BEST EXPRESSION OF GLOBUS DIATHERMY.

TECHNICAL FEATURES

Display

Size 7" color touch screen display

Power

350W $\pm 10\%$

Outputs

2

Emission frequency

400kHz, 470kHz, 700kHz, 1000kHz,
1200kHz, $\pm 10\%$

Treatment modality

Capacitive and resistive

Capacitive electrodes

3 size \varnothing (mm 30, 50, 70)

Resistive electrodes

3 size \varnothing (mm 30, 50, 70)

Mains supply

100-240V AC, 50-60Hz, 1.8A

Size


mm 515x340x170

Weight

gr 10000

EQUIPMENT

- 1 DIACARE 7000 tecartherapy unit
- 1 Counter electrode (plate)
- 1 Monopolar resistive handpiece with 3 heads \varnothing (mm 30, 50, 70)
- 1 Monopolar capacitive handpiece with 3 heads \varnothing (mm 30, 50, 70)
- 1 Bipolar headpiece \varnothing (mm 50)
- 1 Kit for automatic treatment
- 1 Cream can ml 1000
- 1 Power supply unit
- 1 Operating manual



DIACARE 7000

TECARTHERAPY

PRE-SET PROGRAMS TO TREAT THE MOST COMMON PATHOLOGIES

BONE AND ARTICULAR PATHOLOGIES

PATHOLOGIES OF TENDONS AND LIGAMENTS

MUSCULAR AND SOFT TISSUE PATHOLOGIES

VASCULAR AND LYMPHATIC PATHOLOGIES

PATHOLOGIES OF THE NERVOUS SYSTEM



RESISTIVE



CAPACITIVE



BIPOLAR



MULTI-FREQUENCY



PULSED
TREATMENTS



TRANSDERMAL
VEHICULATION



DISPLAY
TOUCH SCREEN



HUMAN ATLAS



HANDS-FREE MODE



PRESET PROGRAMS



AUTOMATIC FREQUENCY
SCANNING



ATHERMAL
TREATMENTS



COLOR DISPLAY



MAINS SUPPLY



DOMINO S.R.L.
Via Vittorio Veneto, 52
31013 Codognè (TV) - ITALY



(+39) 0438 7933



globuscorporation.com

GLOBUS
ITALIAN EXCELLENCE