

# whirlpools



THE STANDARD IN HYDROTHERAPY



A DIVISION OF ACORN ENGINEERING CO.

# quality, function & durability


## THE BEAUTY OF HYDROTHERAPY

Since 1946, Whitehall Manufacturing has provided the health care industry with products unsurpassed in quality, function and durability. Our extensive experience in this industry has given us great insights into the type of quality care necessary for patient rehabilitation. This experience reflects the way we design, manufacture and supply our products.



Our whirlpools' functional features include:






- Powerful 1/2 horsepower turbine motor, UL & CSA approved, is standard on all whirlpool models.
- Stainless steel, (heavy gage, type 304) construction.
- Reinforced stainless steel tank rims.
- Seamless covered bottoms that minimize bacterial build-up.
- Anti-rattle retainer for thermometer.
- Dual scale (Fahrenheit and Celsius) thermometer with adjustable calibration.
- Mobile or stationary designs available.



COVERED BOTTOM



# sizing chart

Model Number/Type	Capacity Gals. (LTRS.)	Tank Dimensions (L x W x D (mm's))	Suggested Use	
P-10-M (Mobile)	10 (38)	22" x 13" x 12" (559 x 330 x 305)		podiatry whirlpools
P-10-MH (Mobile w/Handle)				
P-10-MU (Mobile w/Undercarriage)				
P-10-S (Stationary)				
P-15-M (Mobile)	15 (57)	25" x 13" x 15" (635 x 330 x 381)		podiatry whirlpools
P-15-MH (Mobile w/Handle)				
P-15-MU (Mobile w/Undercarriage)				
P-15-S (Stationary)				
P-22-M (Mobile)	22 (83)	28" x 15" x 18" (711 x 381 x 457)		
P-22-S (Stationary)				
E-15-M (Mobile)	15 (57)	25" x 13" x 15" (635 x 330 x 381)		extremity whirlpools
E-15-MH (Mobile w/Handle)				
E-15-MU (Mobile w/Undercarriage)				
E-15-S (Stationary)				
E-22-M (Mobile)	22 (83)	28" x 15" x 18" (711 x 381 x 457)		extremity whirlpools
E-22-MU (Mobile w/Handle)				
E-22-S (Stationary)				
E-22-SP (Stationary w/Pedestal)				
E-27-M (Mobile)	27 (102)	28" x 15" x 21" (711 x 381 x 533)		
E-27-S (Stationary)				
E-45-M (Mobile)	45 (170)	32" x 15" x 25" (813 x 381 x 635)		
E-45-S (Stationary)				
H-60-M (Mobile)	60 (227)	36" x 20" x 28" (914 x 508 x 711)		hi-boy whirlpools
H-60-S (Stationary)				
H-75-M (Mobile)	75 (284)	42" x 20" x 28" (1067 x 508 x 711)		hi-boy whirlpools
H-75-S (Stationary)				
H-90-M (Mobile)	90 (341)	48" x 20" x 28" (1219 x 508 x 711)		hi-boy whirlpools
H-90-S (Stationary)				
H-105-M (Mobile)	105 (397)	48" x 24" x 28" (1219 x 610 x 711)		
H-105-S (Stationary)				
L-75-M (Mobile)	75 (284)	52" x 24" x 18" (1321 x 610 x 457)		lo-boy whirlpools
L-75-S (Stationary)				
L-90-M (Mobile)	90 (341)	60" x 24" x 18" (1524 x 610 x 457)		lo-boy whirlpools
L-90-S (Stationary)				
L-105-M (Mobile)	105 (397)	66" x 24" x 18" (1676 x 610 x 457)		lo-boy whirlpools
L-105-S (Stationary)				
L-105-SL (Stationary w/Legs)				
S-85-M (Mobile)	85 (322)	48" x 20" x 25" (1219 x 508 x 635)		sports whirlpools
S-85-S (Stationary)				
S-85-SL (Stationary w/Legs)				
S-90-M (Mobile)	90 (341)	46" x 24" x 25" (1168 x 610 x 635)		sports whirlpools
S-90-S (Stationary)				
S-90-SL (Stationary w/Legs)				
S-110-M (Mobile)	110 (416)	56" x 24" x 25" (1422 x 610 x 635)		sports whirlpools
S-110-S (Stationary)				
S-110-SL (Stationary w/Legs)				
SB-100-M (Mobile)	100 (379)	76" x 40" x 21" (1930 x 1016 x 533)		
SB-100-S (Stationary)				
F-260-S (Stationary)	260 (984)	95" x 41" x 22" (2413 x 1041 x 559)		full body immersion tanks
S-300-S (Stationary)				
F-425-S (Stationary)	425 (1,609)	106" x 77" x 22" (2692 x 1956 x 559)		full-body immersion tanks
Hubbard Tank				
S-270-S (Stationary)	270 (1,022)	95" x 41" x 22" (2413 x 1041 x 559)		full-body immersion tanks
Mini-Hubbard Tank				
HYL-10-M (Mobile)	10 (38)	22" x 13" x 12" (559 x 330 x 305)		hydrolift
HYL-15-M (Mobile)	15 (57)	25" x 13" x 15" (635 x 330 x 381)		hydrolift
HYL-22-M (Mobile)	22 (83)	28" x 15" x 18" (711 x 381 x 457)		hydrolift



# podiatry whirlpools

## P series

Podiatry whirlpools provide the ultimate in hydrotherapy treatment for the foot and ankle areas. The powerful 1/2 horsepower turbine motor creates a soothing, cleansing and relaxing treatment for the patient which helps their rehabilitation. The whirlpools are available in several sizes to fit your patient's needs.



P-10-M



P-10-MH (with handle)

# extremity whirlpools

## E series

Versatile and economical best describes the Extremity line of whirlpools. Their compact design saves space and the low gallon capacity reduces hot water usage without sacrificing proper treatment. Each unit treats specific areas of the body such as the foot, ankle, lower leg and knee as well as upper extremities like the elbow, forearm and hand.

The 10, 15, and 22 gallon capacity whirlpools are available with a raised undercarriage and can be utilized for treatment of both upper and lower extremities. The Extremity Series is available in 15, 22, 27 and the larger 45 gallon capacity knee unit, and are available in both mobile and stationary styles.



E-15-M



E-22-MU



E-45-S



# hi-boy whirlpools

## H series

The Hi-Boy Series whirlpools offer complete coverage to the back, hip and shoulder areas. If full immersion is not desired, an optional Adjustable Suspension Seat (-AS) can be utilized. The Hi-Boy Series, when coupled with the optional Tank Top Seat (-TTS), can also be used to treat the lower extremities such as the lower leg, ankle, and foot areas. These whirlpools are available in 60, 75, 90, and 105 gallon capacities in both mobile and stationary styles.



H-75-M



H-105-M



H-105-S

# lo-boy whirlpools

## L series

These long and shallow units are designed for easy access. The lower height makes the Lo-Boy ideal for pediatric and geriatric patients requiring hydrotherapy treatment. The Lo-Boy can provide treatment of various extremity regions along hip and back areas. With the addition of an optional Head Rest (-HR), the patients upper back, shoulder and neck areas can be treated. The Lo-Boy Series is available in 75, 90, and 105 gallon capacities in both mobile and stationary styles.



WHITEHALL MANUFACTURING  
[www.whitehallmfg.com](http://www.whitehallmfg.com)

5



L-75-M-EG



L-105-S

# sports whirlpools

## S series

Based on the needs expressed to us by athletic trainers across the nation, Whitehall has created a full line of whirlpools specifically designed for the treatment and prevention of sports related injuries. These units will facilitate the treatment of the extremities, hips, back, and shoulder areas, plus they are wide enough to accommodate larger patients. The Sports Series is available in 85, 90, and 110 gallon capacities including one style furnished with stainless steel legs (-SL), for easier training room cleaning.

WHITEHALL MANUFACTURING  
www.whitehallmfg.com

6



S-85-S



S-110-SL



SB-100-M



# full body immersion tanks

## F series

Completing our line of whirlpools is the large Full Body Immersion Tank Series. These tanks are designed for burn victim treatment and full body immersion that requires room for body motion and exercise. The two rectangular tanks are available in 260 or 300 gallon capacities. The 425 gallon Hubbard Tank, with its unique figure eight shape, allows therapists easy access to patients. Tanks are stationary in design and come standard with either, twin or single turbine assemblies, a thermometer, head rest and a seamless covered bottom and a complete line of accessories (please reference options chart).



F-260-S



F-425-S



# hydrolift

## HYL series

The perfect apparatus for the raising and lowering of a single whirlpool to treat lower and/or upper extremities. The Hydrolift is fully adjustable, raising the whirlpool rim to any height up to 32", 36" or 38" depending on the model being lifted. In normal operating conditions, the Hydrolift will remain as a stationary lift. However, both Hydrolift and whirlpool may be moved, by lowering the whirlpool to the floor and using the whirlpool casters. A great fit with 10, 15 and 22 gallon capacity whirlpools.



HYL-15-M



HYL-10-M

# enviro-glaze® whirlpools

## -EG option

Whitehall's optional Enviro-Glaze® whirlpool finish can be a colorful and functional alternate to a stainless steel finish. Enviro-Glaze® is a durable, electrostatically sprayed powder, applied to the outside of the tank and then cured at temperatures in excess of 400°F. The finish then becomes an integral part of the stainless steel surface and will not bubble or peel. Our Enviro-Glaze® finish creates a beautiful corrosion and chip-resistant surface that is easily cleaned, maintained and kept in a sanitary condition. Whether you're complimenting an existing décor, matching team colors or a logo, or simply desire a colorful alternative to a stainless steel finish, our colors will give you choices not available anywhere in the industry.





# options




- PHS-2 Patient Handling System with Chair Lift



- OEH Overhead Electric Hoist



- PHS-1 Patient Handling System with Stretcher Lift



-ARS Arm Rest Support



- TTS Tank Top Seat




- CTS Combination Table with Seats



- AS Adjustable Suspension Seat



-HR Headrest, Adjustable



-RS Reclining Seat



-MXT Thermostatic Mixing Valve



HYDRO-CHLOR  
An effective whirlpool  
sanitizer for the treatment  
of burns and open wounds.

# hot & cold therapy overview

## GLACIER FREEZE

A convenient and reliable storage for Glacier Packs. The GlacierFreeze is designed to keep your cold packs properly stored and ready to use.



## THERMALATOR

Whitehall has designed a simple, economical, yet reliable Hot Pack Heating Unit for heating hot packs.

The units are factory set at approximately 168°F, and are constructed completely of stainless steel and come standard with a submersion basket and canvas covered hot packs.



## GLACIER PACKS


The Glacier Packs are ideal for reducing swelling, relieving pain, and relaxing tight joints. The packs are constructed of a durable 14 mil PVC cover. They are available in various sizes.



## PARATHERAPY

Paraffin treatment is a lot easier with our ParaTherapy paraffin treatment unit.

These units are designed to provide heat therapy efficiently, comfortably and is perfect for home or clinic use.



A DIVISION OF ACORN ENGINEERING CO.

15125 Proctor Ave.  
City of Industry, CA 91746

P.O. Box 3527  
City of Industry, CA 91744 U.S.A.  
800-782-7706 626-968-6681  
Fax 626-855-4862  
[www.whitehallmfg.com](http://www.whitehallmfg.com)

L1000910