Ega BedsideCabinet









Features

It is a bedside table with tray to be used next to a hospital bed.

Structure:

All the elements are corrosion protected with the epoxy coating.

Surface finish shock, light, abrasion and scratch resistant, especially in edges and bottom areas. It does not affect the contact with alcohol, solvents, detergents and food products.

Waterproof, fireproof and high resistance HPL (High Pressure Laminated) panels structure. M2 fireproof classification.

Upper cover with perimeter protector made of phenol material HPL and aluminum rod.

Contact surfaces with rounded edges and corners, Sharp edges and cutting surfaces free, which allows an easy cleaning.

Options:

The bedside table has a continuous drawer with sliding guides and anti-fall cap, and a lower cabinet with two doors and inner shelf with bottle hollow.

Removable inner shelf to improve cleaning Access.

Towel rack.

Hollow to be use as tank or bottle holder, or to be use to keep other objects. 5 positions reclining lectern and 65° tilting degree. Enables a 360° rotation for greater bed approach.

4 double pivoting wheels 75 mm. diameter with double rolling rubber band, independent braking in all of them, lint-free and protection disk in all four corners. The two wheels, right under the lectern, are located out of the bedside table main section to increase its stability.

The colour of the front panels can be chosen to mach them with the rest of the room elements.

Reversible table (drawer and cabinet door can be opened from both sides), makes it easier to place the table in both sides of the bed.

Extendable tipping board to approach to the patient and be able to use it as dining table, which encourage comfort.











Ega Bedside Cabinet



Technical information

Maximum/Minimum overbed table height (A / A')	1160 / 950 ± 5 mm.
Overbed table extensibility (B)	135 ± 5 mm.
Overbed table board dimensions (C x D)	620 x 395 ± 5 mm.
Overbed table angles	0° / 60°
Cabinet dimensions (E x F x G)	410 x 400 x 630 ± 5 mm.
Height to the base of the cabinet (H)	105 ± 5 mm.
Total width of the cabinet (without overbed table) (I)	650 ± 5 mm.
Weight	35 kg.
Safety weight load (cabinet / overbed table)	30 / 10 kg.

Colours



