





Option colours



Features

Bedside table with anodised aluminium structure, modular by side panels with laminated coating, attached glue free that provides an innovative aesthetic and extremely light.

- Upper and low textured ABS covers, easy to clean MDF boards with melamine coating and
- "Monoblock" polypropylene drawer, light, easy to clean and maintain.
- 4 double pivoting wheels with brakes with Ø65 and thermoplastic rubber band. This provides it easy to use, safety and quiet operation.
- · Ventilation grid at the back side. Hollow for magazine and inner shelf. Rounded edges and corners.
- Ergonomic zamak handles.
- Easy cleaning.
- Customizable according to client.







Technical information

Structure dimensions:

Width (A)	400 mm. ±10 mm.
Height (B)	800 mm. ±10 mm.
Depth (C)	500 mm. ±10 mm.
Top drawer height (D)	120 mm. ±10 mm.
Hole height (E)	105 mm. ±10 mm.
Bottom drawer height (F)	440 mm. ±10 mm.
Inner shelf height (G)	220 mm. ±10 mm.











dic 2018 Cód: PCO746199_FichaTecnica_MesillaSalduba_EN