



**100 MX**  
**200 MX**  
**300 MX**



LiteGait 100mx

- LiteGait® partial weight bearing gait & balance trainers
- Fall-free Supported Walking
- Use On & Off the Treadmill
- **NEW** Power lift for Improved efficiency



**BUILT-IN POWER LIFT  
FOR EASE OF USE**

- New Battery-Operated Built-in Lift
- Assist Posture, Balance & Weight Bearing
- Differential Control of Right vs. Left Side
- Walk on Treadmill or Over Ground
- Train Gait under Fall Free Conditions
- Frees Therapist to Facilitate Walking
- Wide Range of Patient Sizes & Impairment
- Ideal for Pre-Gait & Upper Body Mobility



LiteGait 300mx



**FREEDOM FOR THERAPISTS INDEPENDENCE FOR PATIENTS**

# 100 MX 200 MX 300 MX

LiteGait family of products have provided safe and effective environment for walking therapy for over a decade. As the original partial weight bearing therapy device, LiteGait has forever altered gait therapy techniques as well as expectations of outcome. As part of the next generation of LiteGait product line, the MX models combine the affordability and versatility of our pediatric manual systems with the ease of use and functionality of our adult LiteGait I systems.

A unique design allows the power lift mechanism to be built directly into the supportive arms of the device. This preserves all of the manual adjustment available in units such as pediatric units while providing 18 inches power lift. You can bring the patient from sit to stand and change the amount of support with a push of a button. The rechargeable battery operated yet lite weight system provides all the functionality you need without the mess of oil pumps or cranks.

Every system comes with everything you need including one harness that fits a large range of patients. The adjustability built into these system allows for growth potential in a single user case; and it allows facilities to easily treat all their patients with the same unit. Consult the specification table to find the right unit for you.

All LiteGait models are backed by one year replacement warranty, a responsive Service department, and an extremely active Education department. When you purchase a LiteGait, you receive the service, the education and the quality you and you clients deserve.



### Specifications\*

	LG-100MX	LG-200MX	LG-300MX
Max Patient weight (lb)	100	200	300
Max Patient Height (in)	61	73	79
Power Lifts (lb/in)	100/18	200/18	300/18
Unit Max/Min Height (in)	65/24	79/30	85/36
Harness fits Girths (in)	17-33	23-45	27-60
Base Inside Width (in)	25/27	27/30/34	30/34
Overall Length (in)	32	42	45

\* approximate and subject to change without notice

**800.332.walk**  
**WWW.LITEGAIT.COM**  
**SALES@LITEGAIT.COM**

### Mobility Research

P.O. Box 3141 • Tempe, Az. 85280 • phone 480.829.1727 • fax 480.829.0737