



LINDA EVO ENT

LINDA EVO ENT is the chair configuration setting designed for Otorhinolaryngology procedures and audiology purposes, which can be equipped with the following specific accessories:

- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- Swivelling left and right armrests
- ELLE Mobile headrest
- ENT footrest



Standard Options

- Swivelling left and right armrests
- Different headrest models available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Different chair painting colors
- Children pillow

ENT Options


- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- ELLE Mobile double-articulated headrest
- COMFORT wide multi-articulated headrest
- Paper roll holder (roll width 45 cm)
- ENT footrest




IR Remote control



ELLE Mobile headrest



COMFORT headrest



Base with rotation



Swivelling armrests



ENT footrest