

LINDA NEXT ENT

LINDA NEXT ENT is the chair configuration setting designed for otorhinolaryngology procedures and audiology purposes, which can be equipped with the following specific accessories:

- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- Swivelling left and right armrests
- LIGHT multi-articulated headrest



Standard Options


- Independent leg section
- Different headrests available according to specific treatment needs
- Memory or hand made soft foam upholsteries (also with handmade contrast stitching)
- Backrest motorized height adjustment
- Backup Battery
- Different chair painting colors
- Bracket for lamp and/or tray connection
- Children pillow
- Leg section foot protection
- Stabilizing plate

ENT Options

- Rotation system around the vertical axis under the chair base (320°)
- IR remote control
- LIGHT multi-articulated headrest
- Swivelling left and right armrests
- COMFORT wide multi-articulated headrest
- Stainless steel lateral tool holder bars
- Adjustable surgical rotating armrests
- Small side instrument tray with chair connection
- Paper roll holder (roll width 45 cm)



COMFORT headrest



IR Remote control



Base with rotation



LIGHT headrest



Swivelling armrests



Paper roll holder