

SERENITY NEXT

DIAGNOSTIC HYSTEROSCOPY | GYNAECOLOGICAL EXAMINATIONS & DIAGNOSTIC | LASER TREATMENTS
MEDICALLY ASSISTED PROCREATION | OBSTETRICS EXAMINATIONS & MONITORING | PROCTOLOGY
SCREENING & COLPOSCOPY MONITORING | ULTRASOUND DIAGNOSTIC | UROLOGY EXAMINATION

Reassuring design, ergonomics and functionality. Lifting capacity up to 250 kg. Easy treating and advanced features. Maximum comfort for both operators and patients.

Standard features

- Complete symmetry for right-handed or left-handed professional
- Electric movement up, down, Trendelenburg, reverse Trendelenburg and backrest recline
- Paper roll holder 60 cm
- Push button control
- Security stops
- Soft start and soft stop movement
- Stainless steel 2,5 L basin
- Wrapping seamless upholsteryes with antimicrobial feature (different colors available)



Options

- 2 auxiliary electric sockets in the base (300W)
- Back-up battery
- Colposcope holder, with anti-crushing safety feature
- Doctor footrest
- Electrically controlled footrests with arc supports
- Electrically heated seat upholstery
- External wheels Ø 50mm
- Fireproof and antistatic fabric available
- Mobile foot control with 3 user profiles including 3 programmable positions for each profile and slow mode program for precise, ultra-slow patient positioning
- Footrests with arc supports and 3-lever adjustment
- Göpel leg holders with 3 levers
- Göpel leg holders with 2 levers and vertical lowering
- Sampling armrest
- Instrument tray
- Integrated extractable leg-rest
- Magnetic headrest cushion
- Memory Foam upholstery
- Movable patient lateral handles
- Quick-release system for emergency backrest manual lowering
- Round basin with valve for urological purposes
- Tool holder lateral bars



Auxiliary sockets



Colposcope holder



Doctor footrest



Electrical footrest



Extractable leg-rest



Mobile foot control



Footrest arc support with 3 levers



Göpel leg holders with 3 levers



Göpel leg holders with 2 levers and vertical adjustment



Sampling armrest



Instrument tray



Magnetic headrest cushion



Memory Foam upholstery



Movable lateral patient handles



Paper roll



Push button control



Quick Release



Stainless steel 2L basin



Tool holder lateral bars



Urology basin

SERENITY NEXT



SERENITY NEXT SETTINGS

SERENITY NEXT GYNE

Standard chair configuration plus:

- Göpel leg holders with 3 levers
- Movable patient lateral handles
- Mobile foot control
- Magnetic headrest cushion
- Push button panel
- 2 L Stainless Steel basin



SERENITY NEXT OBS

Standard chair configuration plus:

- Göpel leg holders with 3 levers
- Movable patient lateral handles
- Integrated extractable leg-rest
- Magnetic headrest cushion
- Push button panel
- 2 L Stainless Steel basin
- Mobile foot control



SERENITY NEXT URO

Standard chair configuration plus:

- Footrest arc support with 3 levers
- Movable patient lateral handles
- Urological basin
- Magnetic headrest cushion
- Wheels
- Push button panel
- Mobile foot control

