ISEN-Mini VNG

IG/IG Pro

Designed for acute care, bedside examination, inpatient consult, and outpatient assessment for differential diagnosis of acute dizziness/vertigo.

IG



IG Pro



Tests

- Spontaneous test
- Gaze test
- **Saccade test (IG Pro)**
- Dix-Hallpike test
- Positional test
- **Caloric test (IG Pro)**
- Head-shake test
- Custom test

Features

- 6.7-inch high-definition touchscreen, 120FPS high image frame rate
- Enlarged patient eye imaging for easy recording of small-scale nystagmus
- Binocular eye image capture, display, storage, and playback. Built-in 256 GB storage space
- Quick emergency recording with a temporary profile
- Data analysis system available for medical record management and exporting of reports and videos
- Wireless data transferring, wireless screen mirroring, wireless printing
- **Real-time nystagmus display. Optional split-screen: display eye image and/or nystagmus tracing (IG Pro)**







Data analysis

ISEN Tech. & Trading Co, Ltd.

302A Bldg.11, No.3 GAO LI ZHANG Rd. HAIDIAN Dist. BEIJING 100095 CHINA

en.isen.com.cn (International) www.isentech.net (USA)

+86-10-62173678 sales@isen.com.cn

>>> Scan to Order Online (USA)

ISEN-Air Caloric Irrigator

Air C

Fast and accurate heating/cooling Real-time video otoscope.

Features

- Real-time video display via a video otoscope, minimize insufficient irrigation
- Durable design. No refilling cooling water needed. Low maintenance
- Voice reminders when reaching preset time/temperature

Air C





Related

ISEN-VNG



ISEN-VNG (Video electronystagmograph), can be used for discriminating between peripheral vertigo and central vertigo, and then diagnosing which side of vestibular system is abnormal in case of peripheral vertigo.

ISEN-vHIT



Independently assess the function of all six semicircular canals.Fast and effective technique for testing the high-frequency VOR pathway. Help distinguishing between central and peripheral vestibular disorders.

ISEN Tech. & Trading Co，Ltd.

302A Bldg.11，No.3 GAO LI ZHANG Rd. HAIDIAN Dist. BEIJING 100095 CHINA

en.isen.com.cn (International) www.isentech.net (USA)

+86-10-62173678 sales@isen.com.cn

>>> Scan to Order Online (USA)